

Product datasheet for RC200625L2

BCAT2 (NM_001190) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: BCAT2 (NM 001190) Human Tagged Lenti ORF Clone

Tag: mGFP Symbol: BCAT2

Synonyms: BCAM; BCATM; BCT2; HVLI; PP18

Mammalian Cell None

Selection:

Vector: pLenti-C-mGFP (PS100071)

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide


The ORF insert of this clone is exactly the same as(RC200625).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF.

ACCN: NM_001190

ORF Size: 1176 bp



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



BCAT2 (NM_001190) Human Tagged Lenti ORF Clone - RC200625L2

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 001190.2</u>

RefSeq Size:1616 bpRefSeq ORF:1179 bp

Locus ID: 587

 UniProt ID:
 O15382

 Cytogenetics:
 19q13.33

Domains: aminotran 4

Protein Families: Druggable Genome

Protein Pathways: Metabolic pathways, Pantothenate and CoA biosynthesis, Valine, leucine and isoleucine

biosynthesis, Valine, leucine and isoleucine degradation

MW: 44.3 kDa


Gene Summary: This gene encodes a branched chain aminotransferase found in mitochondria. The encoded

protein forms a dimer that catalyzes the first step in the production of the branched chain amino acids leucine, isoleucine, and valine. Multiple transcript variants encoding different

isoforms have been found for this gene. [provided by RefSeq, Sep 2009]



Product images:



Circular map for RC200625L2