

## Product datasheet for **RC200652L4V**

### TMED9 (NM\_017510) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	TMED9 (NM_017510) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TMED9
Synonyms:	GMP25; HSGP25L2G; p24a2; p24alpha2; p25
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_017510
ORF Size:	705 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200652).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_017510.4</a>
RefSeq Size:	1402 bp
RefSeq ORF:	708 bp
Locus ID:	54732
UniProt ID:	<a href="#">Q9BVK6</a>
Cytogenetics:	5q35.3
Domains:	EMP24_GP25L
Protein Families:	Transmembrane



[View online »](#)

**MW:** 27.3 kDa

**Gene Summary:** This gene is a member of a family of genes encoding transport proteins located in the endoplasmic reticulum and the Golgi. A similar gene in mouse is the target of microRNA miR-296, which is part of an imprinted cluster. [provided by RefSeq, Jul 2016]