

Product datasheet for **RC200667L3V**

FKBP1B (NM_054033) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	FKBP1B (NM_054033) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FKBP1B
Synonyms:	FKBP1L; FKBP12.6; OTK4; PKBP1L; PPIase
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_054033
ORF Size:	240 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200667).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_054033.1
RefSeq Size:	1016 bp
RefSeq ORF:	243 bp
Locus ID:	2281
UniProt ID:	P68106
Cytogenetics:	2p23.3
Domains:	FKBP
Protein Families:	Druggable Genome



[View online »](#)

MW: 8.8 kDa

Gene Summary: The protein encoded by this gene is a member of the immunophilin protein family, which play a role in immunoregulation and basic cellular processes involving protein folding and trafficking. This encoded protein is a cis-trans prolyl isomerase that binds the immunosuppressants FK506 and rapamycin. It is highly similar to the FK506-binding protein 1A. Its physiological role is thought to be in excitation-contraction coupling in cardiac muscle. There are two alternatively spliced transcript variants of this gene encoding different isoforms. [provided by RefSeq, Jul 2008]