

## Product datasheet for **RC200678L4V**

### **BTF3 (NM\_001207) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	BTF3 (NM_001207) Human Tagged ORF Clone Lentiviral Particle
Symbol:	BTF3
Synonyms:	BETA-NAC; BTF3a; BTF3b; NACB
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001207
ORF Size:	486 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200678).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001207.4</a> , <a href="#">NP_001198.2</a>
RefSeq Size:	1021 bp
RefSeq ORF:	489 bp
Locus ID:	689
UniProt ID:	<a href="#">P20290</a>
Cytogenetics:	5q13.2
Domains:	NAC
Protein Families:	Transcription Factors



[View online »](#)

**MW:** 17.5 kDa

**Gene Summary:** This gene encodes the basic transcription factor 3. This protein forms a stable complex with RNA polymerase IIB and is required for transcriptional initiation. Alternative splicing results in multiple transcript variants encoding different isoforms. This gene has multiple pseudogenes. [provided by RefSeq, Jul 2008]