

Product datasheet for **RC200681L3V**

CHCHD7 (NM_001011671) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CHCHD7 (NM_001011671) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CHCHD7
Synonyms:	COX23
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001011671
ORF Size:	294 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200681).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001011671.2 , NP_001011671.1
RefSeq Size:	1575 bp
RefSeq ORF:	258 bp
Locus ID:	79145
UniProt ID:	Q9BUK0
Cytogenetics:	8q12.1
MW:	11.5 kDa



[View online »](#)