

Product datasheet for **RC200693L1V**

HERPUD1 (NM_001010989) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | HERPUD1 (NM_001010989) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | HERPUD1 |
| Synonyms: | HERP; Mif1; SUP |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-Myc-DDK (PS100064) |
| Tag: | Myc-DDK |
| ACCN: | NM_001010989 |
| ORF Size: | 1173 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC200693). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001010989.1 |
| RefSeq Size: | 2195 bp |
| RefSeq ORF: | 1173 bp |
| Locus ID: | 9709 |
| UniProt ID: | Q15011 |
| Cytogenetics: | 16q13 |
| Protein Families: | Druggable Genome |
| MW: | 43.4 kDa |



[View online »](#)

Gene Summary:

The accumulation of unfolded proteins in the endoplasmic reticulum (ER) triggers the ER stress response. This response includes the inhibition of translation to prevent further accumulation of unfolded proteins, the increased expression of proteins involved in polypeptide folding, known as the unfolded protein response (UPR), and the destruction of misfolded proteins by the ER-associated protein degradation (ERAD) system. This gene may play a role in both UPR and ERAD. Its expression is induced by UPR and it has an ER stress response element in its promoter region while the encoded protein has an N-terminal ubiquitin-like domain which may interact with the ERAD system. This protein has been shown to interact with presenilin proteins and to increase the level of amyloid-beta protein following its overexpression. Alternative splicing of this gene produces multiple transcript variants encoding different isoforms. The full-length nature of all transcript variants has not been determined. [provided by RefSeq, Jan 2013]