

Product datasheet for **RC200694A1V**

Human PCLAF (NM_014736) AAV Particle

Product data:

Product Type:	AAV Particles
Tag:	Myc-DDK
Symbol:	PCLAF
Synonyms:	KIAA0101; L5; NS5ATP9; OEATC; OEATC-1; OEATC1; p15(PAF); p15/PAF; p15PAF; PAF; PAF15
Mammalian Cell	None
Selection:	
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC200694 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s) TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC GCCGCGATCGCC ATGGTGC G GACTAAAGCAGACAGTGTTCAGGCACTTACAGAAAAGTGGTGGCTGCTCGAGCCCCAGAA AGGTGCTTGGTTCTTCCACCTCTGCCACTAATTCGACATCAGTTTCATCGAGGAAAGCTGAAATAAATA TGCAGGAGGGAACCCCGTTTGCCTGCGCCCACTCCCAAGTGGCAAAAAGGAATTGGAGAATTCTTTAGG TTGTCCCTAAAGATTCTGAAAAAGAGAATCAGATTCTGAAGAGGCAGGAAGCAGTGGCTTAGGAAAAG CAAAGAGAAAAGCATGTCCTTTGCAACCTGATCACACAAATGATGAAAAAGAA ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT ACAAGGATGACGACGATAAGGTTTAA
Protein Sequence:	>RC200694 protein sequence Red=Cloning site Green=Tags(s) MVRTKADSVPGTYRKVVAARAPRKVLGSSTSATNSTSVSSRKAENKYAGGNPVCVRPTPKWQKGIGEFFR LSPKDSEKENQIPEEAGSSGLGKAKRKACPLQPDHTNDEKE TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_014736
ORF Size:	333 bp



Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_014736.4
RefSeq Size:	1515 bp
RefSeq ORF:	336 bp
Locus ID:	9768
UniProt ID:	Q15004
Cytogenetics:	15q22.31
MW:	12 kDa