

Product datasheet for **RC200722L3V**

TNFAIP1 (NM_021137) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | TNFAIP1 (NM_021137) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | TNFAIP1 |
| Synonyms: | B12; B61; BTBD34; EDP1; hBACURD2 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_021137 |
| ORF Size: | 948 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC200722). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_021137.3 |
| RefSeq Size: | 3822 bp |
| RefSeq ORF: | 951 bp |
| Locus ID: | 7126 |
| UniProt ID: | Q13829 |
| Cytogenetics: | 17q11.2 |
| Domains: | BTB, K_tetra |
| Protein Families: | Druggable Genome, Ion Channels: Other |



[View online »](#)

MW: 36.2 kDa

Gene Summary: This gene was identified as a gene whose expression can be induced by the tumor necrosis factor alpha (TNF) in umbilical vein endothelial cells. Studies of a similar gene in mouse suggest that the expression of this gene is developmentally regulated in a tissue-specific manner. [provided by RefSeq, Jul 2008]