

Product datasheet for **RC200732L1V**

Fibrillarin (FBL) (NM_001436) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Fibrillarin (FBL) (NM_001436) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Fibrillarin
Synonyms:	FIB; FLRN; Nop1; RNU3IP1
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_001436
ORF Size:	963 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200732).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001436.2
RefSeq Size:	1135 bp
RefSeq ORF:	966 bp
Locus ID:	2091
UniProt ID:	P22087
Cytogenetics:	19q13.2
Domains:	Fibrillarin
Protein Families:	Stem cell - Pluripotency



[View online »](#)

MW: 33.6 kDa

Gene Summary: This gene product is a component of a nucleolar small nuclear ribonucleoprotein (snRNP) particle thought to participate in the first step in processing preribosomal RNA. It is associated with the U3, U8, and U13 small nuclear RNAs and is located in the dense fibrillar component (DFC) of the nucleolus. The encoded protein contains an N-terminal repetitive domain that is rich in glycine and arginine residues, like fibrillarins in other species. Its central region resembles an RNA-binding domain and contains an RNP consensus sequence. Antisera from approximately 8% of humans with the autoimmune disease scleroderma recognize fibrillarlin. [provided by RefSeq, Jul 2008]