

Product datasheet for **RC200741L4V**

SURF6 (NM_006753) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SURF6 (NM_006753) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SURF6
Synonyms:	RRP14
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_006753
ORF Size:	1083 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200741).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006753.3
RefSeq Size:	2349 bp
RefSeq ORF:	1086 bp
Locus ID:	6838
UniProt ID:	O75683
Cytogenetics:	9q34.2
Domains:	SURF6
MW:	41.3 kDa



[View online »](#)

Gene Summary:

This gene encodes a conserved protein that is localized to the nucleolus. The encoded protein may function as a nucleolar-matrix protein with nucleic acid-binding properties. There is a pseudogene for this gene on chromosome Y. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2013]