

Product datasheet for **RC200745L3V**

KCTD14 (NM_023930) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	KCTD14 (NM_023930) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KCTD14
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_023930
ORF Size:	675 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200745).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_023930.3 , NP_076419.1
RefSeq Size:	1688 bp
RefSeq ORF:	768 bp
Locus ID:	65987
UniProt ID:	Q9BQ13
Cytogenetics:	11q14.1
Domains:	BTB, K_tetra
Protein Families:	Ion Channels: Other
MW:	26.2 kDa



[View online »](#)