

## Product datasheet for **RC200749L3V**

### **MCG10 (PCBP4) (NM\_033008) Human Tagged ORF Clone Lentiviral Particle**

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | MCG10 (PCBP4) (NM_033008) Human Tagged ORF Clone Lentiviral Particle   |
| Symbol:                   | MCG10  |
| Synonyms:                 | CBP; LIP4; MCG10   |
| Mammalian Cell Selection: | Puromycin  |
| Vector:                   | pLenti-C-Myc-DDK-P2A-Puro (PS100092)   |
| Tag:                      | Myc-DDK  |
| ACCN:                     | NM_033008  |
| ORF Size:                 | 1209 bp  |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC200749).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_033008.1</a>  |
| RefSeq Size:              | 2040 bp  |
| RefSeq ORF:               | 1212 bp  |
| Locus ID:                 | 57060  |
| UniProt ID:               | <a href="#">P57723</a>   |
| Cytogenetics:             | 3p21.2   |
| Domains:                  | KH   |
| MW:                       | 41.5 kDa   |


[View online »](#)

**Gene Summary:**

This gene encodes a member of the KH-domain protein subfamily. Proteins of this subfamily, also referred to as alpha-CPs, bind to RNA with a specificity for C-rich pyrimidine regions. Alpha-CPs play important roles in post-transcriptional activities and have different cellular distributions. This gene is induced by the p53 tumor suppressor, and the encoded protein can suppress cell proliferation by inducing apoptosis and cell cycle arrest in G(2)-M. This gene's protein is found in the cytoplasm, yet it lacks the nuclear localization signals found in other subfamily members. Multiple alternatively spliced transcript variants have been described, but the full-length nature for only some has been determined. [provided by RefSeq, Jul 2008]