

Product datasheet for **RC200766A1V**

Human RIBC2 (NM_015653) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human RIBC2 (NM_015653) AAV Particle
Tag:	Myc-DDK
Symbol:	RIBC2
Synonyms:	C22orf11; TRIB
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC200766 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAGGCAAAATGACAAAATCATGTGCATATTGGAAAACCGAAAAAGAGGGATAGGAAAAATCTCTGTA
GGGCTATCAATGACTTCCAACAGAGCTTTCAGAAGCCAGAACTCGCCGTGAATTTGATCTGTCCGACCC
CCTAGCCCTTAAGAAAGATCTTCCAGCCCGGCAGTCAGATAATGATGTTTCGGAATACGATATCAGGAATG
CAGAAATTCATGGGAGAGGATTTAACTTCCATGAGAGGAAGAAATCCAAGAGGAACAAAACAGAGAAT
GGTCTTTGCAGCAGCAAAGGGAATGGAAGAACGCCCGTGTCTGAACAAAAATGCGCAGAGGCCCTCTACAC
AGAGACAAGGCTGCAGTTTGACGAGACAGCCAAGCACCTCCAGAAGCTGGAAAGCACCACAGAAAGGCA
GTTTGTGCATCTGTGAAAGACTTCAACAAGAGCCAGGCCATCGAGTCAGTGAAAGGAAAAAGCAAGAGA
AAAAGCAAGAACAAGAGGACAATTGGCCGAGATCACCAACCTCCTGCGTGGGGACCTGCTCTCCGAGAA
CCCGCAGCAGGCAGCCAGCTCCTTGGGCCCCACCGCTGGTCCCTGACCGCTGGAAGGGCATGACCCAG
GAGCAGCTGGAGCAGATCCGCCTAGTCCAGAAGCAGCAAATCCAGGAGAAGCTGAGGCTCCAGGAAGAAA
AGCGCCAGCGAGACCTGGACTGGGACCGCGGAGGATTCAGGGGGCTCGCGCCACCCTGCTGTTTGAAGCG
GCAGCAGTGGCGGCGCAGCGACCTGCGCAGAGCTCTGGACAGCAGCAACCTCAGCCTGGCCAAGGAG
CAGCATTTGCAGAAAAATATATGAATGAAGTCTATACAAATCAACCCACGGGAGACTATTTACACAAT
TTAATACAGGAAGTCGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC200766 protein sequence Red=Cloning site Green=Tags(s) MRQNDKIMCILENRKKRDRKNLCRAINDFQQSFQKPTRREFDLSDPLALKKDLPARQSDNDVRNTISGM QKFMGEDLNHHERKKFQEEQNREWSLQQREWKNARAEQKCAEALYTETRLQFDETAHLQKLESTTRKA VCASVKDFNKSQAIESVERKKQEKKEQEDNLAEITNLLRGDLLSENPPQAASSFGPHRVVPDRWKGMTQ EQLEQIRLVQKQIQEKLRLQEEKRQRDLWDRRRIQGARATLLFERQQWRRQRDLRRALDSSNLSLAKE QHLQKKYMNEVYTNQPTGDYFTQFNTGSR TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_015653
ORF Size:	927 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_015653.2 , NP_056468.1
RefSeq Size:	1405 bp
RefSeq ORF:	1149 bp
Locus ID:	26150
UniProt ID:	Q9H4K1
Cytogenetics:	22q13.31
MW:	37.1 kDa