

Product datasheet for **RC200767A1V**

Human TSPAN15 (NM_012339) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human TSPAN15 (NM_012339) AAV Particle
Tag:	Myc-DDK
Symbol:	TSPAN15
Synonyms:	2700063A19Rik; NET-7; NET7; TM4SF15
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC200767 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCCGCGCGGGGACTCGGAGCAGGTGCGCTACTGCGCGCGCTTCTCCTACCTCTGGCTCAAGTTTTAC
TTATCATCTATTCCACCGTGTCTGGCTGATTGGGGCCCTGGTCCTGTCTGTGGGCATCTATGCAGAGGT
TGAGCGGCAGAAATATAAAACCTTGAAAGTGCCTTCTGGCTCCAGCCATCATCCTCATCCTCCTGGGC
GTCGTCATGTTTCATGGTCTCCTTCATTGGTGTGCTGGCGTCCCTCCGTGACAACCTGTACCTTCTCCAAG
CATTTCATGTACATCCTTGGGATCTGCCATCATGAGCTCATTGGTGGCGTGGTGGCCTTGACCTTCCG
GAACCAGACCATTGACTTCCTGAACGACAACATTGGAAGAGGAATTGAGAACTACTATGATGATCTGGAC
TTCAAAAACATCATGGACTTTGTTCAAAAAAGTTCAAGTGCTGTGGCGGGGAGGACTACCGAGATTGGA
GCAAGAATCAGTACCACGACTGCAGTGCCCTGGACCCCTGGCCTGTGGGGTGCCTACACCTGCTGCAT
CAGGAACACGACAGAAGTTGTCAACACCATGTGTGGCTACAAAATATCGACAAGGAGCGTTTCAGTGTG
CAGGATGTCATCTACGTGCGGGGCTGCACCAACGCCGTGATCATCTGGTTCATGGACAACACCATCA
TGGCGGGCATCCTCCTGGGCATCCTGCTCCCCAGTTCCTGGGGGTGCTGCTGACGCTGCTGTACATCAC
CCGGGTGGAGGACATCATGAGCACTCTGTCACTGATGGGCTCCTTGGGCCCGGTGCCAAGCCCAGC
GTGGAGGCGGCAGGCACGGGATGCTGCTGTGTGCTACCCCAAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC200767 protein sequence Red=Cloning site Green=Tags(s) MPRGDSEQVRYCARFSYLWKFSLIIYSTVFWLIGALVLSVGIIAEVERQKYKTLESAFLAPAIILILLG VVMFMVSFIGVLASLRDNL YLLQAFMYILGICLIMELIGGVVALTFRNQTIDFLNDNIRRGIIENYYDDL FKNIMDFVQKKFKCCGGEDYRDWSKNQYHDCSAPGPLACGVPTCCIRNTTEVVNTMCGYKTIDKERFSV QDVIYVRGCTNAV IIFWMDNYTIMAGILLGILLPQFLGVLLTLLYITRVEDIIMEHSVTDGLLGPGAKPS VEAAGTGCCLCYPN TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_012339
ORF Size:	882 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_012339.3</u>
RefSeq Size:	1726 bp
RefSeq ORF:	885 bp
Locus ID:	23555
UniProt ID:	<u>O95858</u>
Cytogenetics:	10q22.1
MW:	33.2 kDa