

Product datasheet for **RC200774A1V**

Human RPS15A (NM_001019) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human RPS15A (NM_001019) AAV Particle
Tag: Myc-DDK
Symbol: RPS15A
Synonyms: DBA20; S15a
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC200774 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGTGCGCATGAATGTCCTGGCAGATGCTCTCAAGAGTATCAACAATGCCGAAAAGAGAGGCAAACGCC
 AGGTGCTTATTAGGCCGTGCTCCAAAGTCATCGTCCGGTTTCTCACTGTGATGATGAAGCATGGTTACAT
 TGGCGAATTTGAAATCATTGATGACCACAGAGCTGGGAAAATTGTTGTGAACCTCACAGGCAGGCTAAAC
 AAGTGTGGGTGATCAGCCCCAGATCTGACGTGCAACTCAAAGACCTGGAAAAATGGCAGAATAATCTGC
 TTCCATCCCAGTTTGGTTTCATTGTACTGACAACCTCAGCAGTGGCATCATGGACCATGAAGAAGCAAG
 ACGAAAACACACAGGAGGGAAAATCCTGGGATTCCTTTTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC200774 protein sequence
 Red=Cloning site Green=Tags(s)

MVRMNVLADALKSINNAEKRGKRQVLIRPCSKVIVRFLTVMMKHGYIGEFEIIDDHRAGKIVVNLTGRLN
 KCGVISPRSDVQLKDLEKWQNNLLPSRQFGFIVLTTTSAGIMDHEEARRKHGGKILGFFF

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_001019



[View online »](#)

ORF Size:	390 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001019.4</u>
RefSeq Size:	554 bp
RefSeq ORF:	393 bp
Locus ID:	6210
UniProt ID:	<u>P62244</u>
Cytogenetics:	16p12.3
MW:	14.8 kDa