

Product datasheet for **RC200777L4V**

hnRNP G (RBMX) (NM_002139) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	hnRNP G (RBMX) (NM_002139) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RBMX
Synonyms:	hnRNP-G; HNRNPG; HNRPG; MRXS11; RBMXP1; RBMXRT; RNMX
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_002139
ORF Size:	1173 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200777).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002139.2
RefSeq Size:	2097 bp
RefSeq ORF:	1176 bp
Locus ID:	27316
UniProt ID:	P38159
Cytogenetics:	Xq26.3
Domains:	RRM, RRM_1
Protein Pathways:	Spliceosome



[View online »](#)

MW: 42.3 kDa

Gene Summary: This gene belongs to the RBMY gene family which includes candidate Y chromosome spermatogenesis genes. This gene, an active X chromosome homolog of the Y chromosome RBMY gene, is widely expressed whereas the RBMY gene evolved a male-specific function in spermatogenesis. Pseudogenes of this gene, found on chromosomes 1, 4, 9, 11, and 6, were likely derived by retrotransposition from the original gene. Alternatively spliced transcript variants encoding different isoforms have been identified. A snoRNA gene (SNORD61) is found in one of its introns. [provided by RefSeq, Sep 2009]