

Product datasheet for **RC200780A1V**

Human ZNF747 (NM_023931) AAV Particle

Product data:

Product Type: AAV Particles
Tag: Myc-DDK
Symbol: ZNF747
Mammalian Cell None
Selection:
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC200780 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGACCGATCCGTCGCTGGGGCTGACAGTCCCCATGGCGCCGCTCTGGCCCCGCTCCCTCCCCGGGACC
CAAACGGGGCGGGATCCGAGTGGAGAAAGCCCGGGCCGTGAGCTTCGCGACGTGGCCGTGACTTCTC
CCGGGAGGAGTGGGGTGCCTGCGGCCCGCGCAGAGGGCCCTGTACCGGACGTGATCGGGAGACCTAC
GGCCACCTGGGCGCGCTCGGTGAGAGCCCACTGCTTGCTGGGCCCTGCGCCTCCACAGGCCCTGCCG
CGCCTCTGGGAGCTGCGTGTGGAGTTGGGGGCCCGGGGCCGGGACGGCGGCCTCCTCGCAGCGTGGGGT
TTGCGTTCTTCTCCCCAGGAGTCGGAGGCAGCAAGCCGGCGCTCATCTCTGGGTGGAGGAGAAGGCCG
AACTGTGGGATCCGGCTGCCAGGATCCGGAGGTGGCGAAGTGTCCGACAGAAGCGGACCCAGCAGATTC
CAGAAACAAGGAAGGAAAGACAAAGGAAGGGACGGGAGCCCTGGAGAAGCCCGACCTGTGGCCGCC
GGGTCTCCTGGGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC200780 protein sequence
Red=Cloning site Green=Tags(s)

MTDPSLGLTVPMAPPLAPLPPRDPNGAGSEWRKPGAVSFADVAVYFSREEWGCLRPAQRALYRDVMRETY
GHLGALGESPTCLPGPCASTGPAAPLGAACGVGGPGAGQAASSQRGVCVLLLPQESEAASRRSSPGWRRRP
NCGIRLPRIIRWRSVRQKRTQQIPETRRKDKGKGREPWSPTLWPPGLLG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2



ACCN:	NM_023931
ORF Size:	573 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_023931.2</u>
RefSeq Size:	2392 bp
RefSeq ORF:	576 bp
Locus ID:	65988
UniProt ID:	<u>Q9BV97</u>
Cytogenetics:	16p11.2
MW:	20.6 kDa