

Product datasheet for **RC200782A1V**

Human C11orf48 (LBHD1) (NM_024099) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human C11orf48 (LBHD1) (NM_024099) AAV Particle
Tag:	Myc-DDK
Symbol:	C11orf48
Synonyms:	C11orf48
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC200782 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCCCTTGTGCCAGGGAGAAGCAAGGAGGATGGGCTTTGGACTAGAAATAGCCCAGGCTCCTCCCAGC
ATCCAGAAAGTCCCAGGCTGCCCAACCTCTCTGGGACAGAGGAAAAATTGGCAAGGTTGAAGGTCACCA
GCACATTCAGGATTTCTCTCAAAGTCCCATCTGCCGTCTATTGGTGGATCCAGTGAGGTGAATGAA
GAGAGTGGGGATCTCCATTTGCCCCATGAGGAGCTGCTGCTGCTCACTGATGGTGAGGAAGAGGATGCTG
AGGCCTTCTTCAAGACCAAAGTGAAGAGCCAGGCTGGGCTTGGAGCCACAGGACCCTAGAAGTCCTTT
AAGAACATTTAACGCTGGACTCAGCTGGGGGCAGGACCAGGATGAAGAAGATGCTTGTGGATTCTTGAG
GACACAGCATGTCTGGAAGCCACCAACCACTGTCCCTTCTGGGACTCAACAGGCTCCCGTGTGTGTAGAA
GTGGCTTTGTGGAATATTCCTATCTCTGCCTCCTAATAGCTTTGAGGGAGCTGAAGAAGAAGCTGTTCA
AACGCCGGCGGGTGTGAATCGGGAGCGGCGTCTGAGGCACCGGGTGGTGGGGCTGTGATAGACCAAGG
GCTGATCACGGGCACCACTCAAGAAGCGGGCGTCCAGTGCACGTGCCAACATTACACTGTCAGGGAAG
AAGCGCAGAAACTCTCCAGCAGATCCGGCTTGCCAGAAAGAGAAGACAGCCATGGAAGTGAAGCCC
CTTCAAAGCCAGCCAGGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC200782 protein sequence Red=Cloning site Green=Tags(s) MALVPGRSKEDGLWTRNSPGSSQHPESPRLPNPLWDRGKIGKVEGHQHIQDFSQKSHLPSIVVESSEVNE ESGDLHLPHEELLLLTDGEEEDAEAFFQDQSEEPGWAWSPQDPRSPLRTFNAGLSWGQDQDEEDACWILE DTACLEATNHCPFWSTGSRVCRSGFVEYSHLLPPNSFEGAEAAVQTPAGVESGAASEAPGGRGCDRPR ADHAAPPQEAGVQCTCQHYTVREEAQKTPADPACPEREDSHGSGSPFKASQD TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_024099
ORF Size:	789 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_024099.3, NP_077004.2</u>
RefSeq Size:	1357 bp
RefSeq ORF:	792 bp
Locus ID:	79081
UniProt ID:	<u>Q9BQE6</u>
Cytogenetics:	11q12.3
MW:	28.8 kDa