

Product datasheet for **RC200783**

PPDPF (NM_024299) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids

Tag: Myc-DDK

Symbol: PPDPF

Synonyms: C20orf149; dJ697K14.9; exdpf

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >RC200783 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGCGGCCATCCCTCCAGCGGCTCGCTCGTGCCACCCACGACTACTACCGGCGCCGCTGGGTTCCA
CTTCCAGCAACAGCTCCTGCAGCAGTACCGAGTGCCCCGGGAAGCCATTCCCCACCCCCAGGTCTCCC
CAAGGCTGACCCGGGTCAATGGTGGGCCAGCTTCTTTTCGGGAAGTCCACCTCCCGTTTCATGGCCACG
GTGTTGGAGTCCGAGAGCACTCGGAACCTCCCAGGCCTCCAGCAGCATGACCGCTGTGGCTGGCTC
GGGACGCCCGGAGGAAGCAGCCGCGGTCACTCCAGCACAGCCAGCGCTGGGCCCCCGTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC200783 protein sequence
Red=Cloning site Green=Tags(s)

MAAIPSSGSLVATHDYRRRLGSTSSNSSCSSTECPEAIPHPPLPKADPGHWWASFFFGKSTLPFMT
VLESAEHSEPPQASSMTACGLARDAPRKQPGGQSSTASAGPPS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Chromatograms: https://cdn.origene.com/chromatograms/mk6386_b09.zip

Restriction Sites: SgfI-MluI



Cloning Scheme:

Cloning sites used for ORF Shuttling:



ACCN: NM_024299

ORF Size: 342 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10uq of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22µm filter is required.

RefSeq: NM_024299.4

RefSeq Size: 794 bp

RefSeq ORF: 345 bp

Locus ID: 79144


UniProt ID: [Q9H3Y8](#)

Cytogenetics: 20q13.33


MW: 11.8 kDa

Gene Summary: Probable regulator of exocrine pancreas development.[UniProtKB/Swiss-Prot Function]


Product images:



Circular map for RC200783



Western blot validation of overexpression lysate (Cat# [LY411320]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from un-transfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC200783 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified PPDPF protein (Cat# [TP300783]). The protein was produced from HEK293T cells transfected with PPDPF cDNA clone (Cat# RC200783) using MegaTran 2.0 (Cat# [TT210002]).