

## Product datasheet for **RC200786L3V**

### POFUT1 (NM\_172236) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	POFUT1 (NM_172236) Human Tagged ORF Clone Lentiviral Particle
Symbol:	POFUT1
Synonyms:	DDD2; FUT12; O-Fuc-T; O-FucT-1; O-FUT; OFUCT1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_172236
ORF Size:	582 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200786).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_172236.1</a>
RefSeq Size:	1524 bp
RefSeq ORF:	585 bp
Locus ID:	23509
UniProt ID:	<a href="#">Q9H488</a>
Cytogenetics:	20q11.21
Protein Families:	Druggable Genome, Transcription Factors
MW:	22.3 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a member of the glycosyltransferase O-Fuc family. This enzyme adds O-fucose through an O-glycosidic linkage to conserved serine or threonine residues in the epidermal growth factor-like repeats of a number of cell surface and secreted proteins. O-fucose glycans are involved in ligand-induced receptor signaling. Alternative splicing of this gene results in two transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]