

## Product datasheet for **RC200797**

### ZNF576 (NM\_024327) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** ZNF576 (NM\_024327) Human Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** ZNF576  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >RC200797 ORF sequence  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGAGGACCCGAACCTGAAGAGAACATGAAGCAGCAGGATTCACCCAAGGAGAGAAGTCCCCAGAGCC  
CAGGAGGCAACATCTGCCACCTGGGGGCCCGAAGTGCACCCGCTGCCTCATCACCTTCGAGATTCCAA  
GTTCCAGGAGCGTCACATGAAGCGGGAGCACCCAGCGGACTTCGTGGCCAGAAGCTGCAGGGGGTCCCT  
TTCATCTGCTTACCTGCGCCCGCTCCTTCCTCCTCCAAAGCCCTAATCACCCACCAGCGCAGCCAGC  
GTCCAGCCGCAAGCCACCTGCCGTTGCAACCACTACTGCCAGCCACCTCCCTTGTCTGCTGACTG  
TGGCAAGACCTTTGGGCAGGCTGTTTCTCTGAGGCGGCACCGCCAGATGCATGAGGTCCGTGCCCTCCT  
GGCACCTTCGCTGCACAGAGTGCAGTTCAGGACTTTGCTCAGGAAGCAGGCTGCATCAACACTACATTC  
GGCATGCCCGGGGGAGCTC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC200797 protein sequence  
 Red=Cloning site Green=Tags(s)

MEDPNPEENMKQQDSPKERSPQSPGGNICHLGAPKCTRCLITFADSKFQERHMKREHPADFVAQKLQGVL  
FICFTCARSLSSKALITHQRSHGPAAKPTLPVATTTAQTFFPCPDCGKTFGQAVSLRRHRQMHEVRAPP  
GTFACTECGQDFAEQEAGLHQHYIRHARGEL

**TRTRPLEQKLI**SEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mk6394\\_c12.zip](https://cdn.origene.com/chromatograms/mk6394_c12.zip)

**Restriction Sites:** SgfI-MluI



View online »

**Cloning Scheme:**


**ACCN:** NM\_024327

**ORF Size:** 510 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_024327.2](#)

**RefSeq Size:** 1899 bp

**RefSeq ORF:** 513 bp

**Locus ID:** 79177

**UniProt ID:** [Q9H609](#)


**Cytogenetics:** 19q13.31

Domains: zf-C2H2


MW: 18.9 kDa

Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for RC200797



Western blot validation of overexpression lysate (Cat# [LY428834]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with [RC227781] using transfection reagent MegaTran 2.0 (Cat# [TT210002]).