

## Product datasheet for **RC200801L4V**

### SPSB2 (NM\_032641) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SPSB2 (NM_032641) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SPSB2
Synonyms:	GRCC9; SSB2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_032641
ORF Size:	789 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200801).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_032641.1</a>
RefSeq Size:	1272 bp
RefSeq ORF:	792 bp
Locus ID:	84727
UniProt ID:	<a href="#">Q99619</a>
Cytogenetics:	12p13.31
Domains:	SOCS, SPRY
Protein Families:	Druggable Genome



[View online »](#)

**MW:** 28.6 kDa

**Gene Summary:** This gene encodes a member of a subfamily of proteins containing a central SPRY (repeats in splA and RyR) domain and a C-terminal suppressor of cytokine signaling (SOCS) box. This protein plays a role in cell signaling. This gene is present in a gene-rich cluster on chromosome 12p13 in the vicinity of the CD4 antigen and triosephosphate isomerase genes. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2016]