

Product datasheet for **RC200805A1V**

Human KLF13 (NM_015995) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human KLF13 (NM_015995) AAV Particle
Tag:	Myc-DDK
Symbol:	KLF13
Synonyms:	BTEB3; FKLF2; NSLP1; RFLAT-1; RFLAT1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC200805 representing NM_015995 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCAGCCGCCCTATGTGGACCACTTCGCCGCCGAGTGCCTCGTGTCCATGTCGAGCCGCGCGGTGC
TGCACGGGCGCGGGAGGGGCCGAGTCCCGGCCGAGGGCGCGGCGGTGGCCGCCACCCCCACGCTGCC
CCGCGTCGAGGAGCGCCGCGACGTAAGGACAGCGCCTCGCTCTTCGTGGTGGCGCGGATCCTAGCGGAC
CTCAACCAGCAAGCGCCGCGCCCGCCCGCGGAGCGCAGGGAGGGCGCGCGGCCCGGAAGGCGAGGA
CCCCCTGCCGCTGCCGCGCCCGCCCGAGCCACCTCCCCGCGGCCGAAGGCGCGCGCGCCGCGCC
CCCCAGCCCGCGGTGGAGCGAGCCGAGCCCGAGGCGGGGCTGGAGCCCGAGCGGGAGCCGGGGCCGCG
GGGAGCGGCGAGCCCGCCTCAGACAAAGGTCGCGCGGGGCCGAAGTCGCGCCGACCTCGAGTCCCCGC
AGAGGAAGCACAAAGTGCCACTACGCGGGCTGCGAGAAAGTTTACGGGAAATCTTCGCACCTCAAGGCGCA
CCTGAGAACTCACACAGGTGAGAGGCCCTTCGCTGCAGCTGGCAGGACTGCAACAAGAAGTTCGCGCGC
TCCGACGAGCTGGCGCGCACTACCGCACACACGCGGCGAGAAGAAGTTCAGCTGCCCCATCTGCGAGA
AGCGCTTCATGCGCAGCGACCACTGACCAAGCACGCGCGCCGACGCCAACTTCCACCCGGGAATGCT
GCAGCGGCGCGCGGGGCTCGCGGACCGGCTCCCTCAGCGACTACAGCCGCTCCGACGCCAGCAGCCCC
ACCATCAGCCCGCCAGCTCGCCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC200805 representing NM_015995 Red=Cloning site Green=Tags(s) MAAAAYVDHFAAECLVSMSSRAVVHGPREGPESRPEGAAVAATPTLPRVEERRDGKDSASLFVVARILAD LNQQAPAPAPAAERREGAAARKARTPCRLPPAPEPTSPGAEGAAAAPPSPAWSEPEPEAGLEPEREPGPA GSGEPGLRQRRGRSRADLESPQRKHKCHYAGCEKVYKSSHLKAHLRTHTGERPFACSWQDCNKKFAR SDELARHYRTHTGEKKFSCPICEKRFMRSDHLTKHARRHANFHPGMLQRRGGGSRTGSLSDYSRSDASSP TISPASSP TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_015995
ORF Size:	864 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_015995.2</u>
RefSeq Size:	6817 bp
RefSeq ORF:	867 bp
Locus ID:	51621
UniProt ID:	<u>Q9Y2Y9</u>
Cytogenetics:	15q13.3
MW:	31 kDa