

Product datasheet for **RC200811A1V**

Human FSTL3 (NM_005860) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human FSTL3 (NM_005860) AAV Particle
Tag:	Myc-DDK
Symbol:	FSTL3
Synonyms:	FLRG; FSRP
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC200811 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCGTCCCGGGCGCCAGGGCCACTCTGGCCTCTGCCCTGGGGGGCCCTGGCTTGGGCCGTGGGCTTCG
TGAGCTCCATGGGCTCGGGGAACCCCGCGCCCGGTGGTGTGCTGGCTCCAGCAGGGCCAGGAGGCCAC
CTGCAGCCTGGTGCTCCAGACTGATGTACCCGGGCGAGTGCTGTGCCTCCGGCAACATTGACACCGCC
TGGTCCAACCTCACCCACCGGGGAACAAGATCAACCTCCTCGGCTTCTTGGGCTTGTCCACTGCCTTC
CCTGCAAGATTCTGTCGACGGCGTGGAGTGCGGCCCGGGCAAGGCGTGCCGCATGCTGGGGGGCCGCC
GCGCTGCGAGTGCGCGCCGACTGCTCGGGGCTCCCGGCGCGGTTGCAGGTCTGCGGCTCAGACGGCGCC
ACCTACCGCGACGAGTGCGAGTGCGCGCCGCGCGCTGCCGCGGCCACCCGACCTGAGCGTCATGTACC
GGGGCCGCTGCCGCAAGTCTGTGAGCACGTGGTGTGCCCGCGGCCACAGTCGTGCGTCGTGGACCAGAC
GGGCAGCGCCCACTGCGTGGTGTGTCGAGCGCGCCCTGCCCTGTGCCCTCCAGCCCCGCCAGGAGCTT
TGCGGCAACAACAGTCACCTACATCTCCTGTGCCACATGCCCAGGCCACCTGCTTCTGGGCGCT
CCATCGGCGTGCGCCACGCGGGCAGTGCGCAGGCACCCCTGAGGAGCCGCCAGGTGGTGTGAGTCTGCAGA
AGAGGAAGAGAACTTCGTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC200811 protein sequence Red=Cloning site Green=Tags(s) MRPGAPGPLWPLPWGALAWAVGFVSSMSGNPAPGGVCWLQGGQEATCSLVLQTDVTRAECASGNIDTA WSNLTHPGNKINLLGFLGLVHCLPCKDSCDGVCEGPGKACRMLGGRPRCECAPDCSGLPARLQVCGSDGA TYRDECELRAARCRGHPDLVSMYRGRCRKSCHEVVCPRPQSCVVDQTGSAHCVVCRAAPCPVPSSPGQEL CGNNNVTYISSCHMRQATCFLGRSIGVRHAGSCAGTPEEPPGGESAEEEEENFV TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_005860
ORF Size:	789 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_005860.2</u>
RefSeq Size:	2525 bp
RefSeq ORF:	792 bp
Locus ID:	10272
UniProt ID:	<u>O95633</u>
Cytogenetics:	19p13.3
MW:	27.7 kDa