

## Product datasheet for **RC200811L3V**

### **FSTL3 (NM\_005860) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	FSTL3 (NM_005860) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FSTL3
Synonyms:	FLRG; FSRP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_005860
ORF Size:	789 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200811).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_005860.2</a>
RefSeq Size:	2525 bp
RefSeq ORF:	792 bp
Locus ID:	10272
UniProt ID:	<a href="#">O95633</a>
Cytogenetics:	19p13.3
Domains:	kazal, FOLN
Protein Families:	Secreted Protein


[View online »](#)

**MW:** 27.7 kDa

**Gene Summary:** Follistatin-like 3 is a secreted glycoprotein of the follistatin-module-protein family. It may have a role in leukemogenesis. [provided by RefSeq, Jul 2008]