

Product datasheet for **RC200819A1V**

Human FGFBP1 (NM_005130) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human FGFBP1 (NM_005130) AAV Particle
Tag:	Myc-DDK
Symbol:	FGFBP1
Synonyms:	FGF-BP; FGF-BP1; FGFBP; FGFBP-1; HBP17
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC200819 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAAGATCTGTAGCCTCACCTGCTCTCCTTCTCTACTGGCTGCTCAGGTGCTCCTGGTGGAGGGGA
AAAAAAAGTGAAGAATGGACTTCACAGCAAAGTGGTCTCAGAACAAAGGACACTCTGGGCAACACCCA
GATTAAGCAGAAAAGCAGGCCCGGGAACAAAGGCAAGTTTGTACCAAAGACCAAGCCAAGTGCAGATGG
GCTGCTACTGAGCAGGAGGAGGGCATCTCTCAAGTTGAGTGCACTCAATTGGACCATGAATTTTCCT
GTGCTTTTGTGGCAATCCAACCTCATGCCTAAAGCTCAAGGATGAGAGAGTCTATTGGAACAAGTTGC
CCGGAATCTGCGCTCACAGAAAGACATCTGTAGATATTCCAAGACAGCTGTGAAAACAGAGTGTGCAGA
AAGGATTTTCCAGAATCCAGTCTTAAGCTAGTCAGCTCCTCTATTTGGGAACACAAAGCCCAGGAAGG
AGAAAACAGAGATGTCCCCAGGGAGCACATCAAGGGCAAAGAGACCACCCCTCTAGCCTAGCAGTGAC
CCAGACCATGGCCACCAAGCTCCCGAGTGTGTGGAGGACCCAGATATGGCAAACCAGAGGAAGACTGCC
CTGGAGTTCTGTGGAGAGACTTGGAGCTCTCTGCACATTCTCCTCAGCATAGTGCAGGACACGTCAT
GC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC200819 protein sequence Red=Cloning site Green=Tags(s) MKICSLTLLSFLLLAAQVLLVEGKKVKNGLSKVVSEQKDTLGNTQIKQKSRPGNKGKFVTKDQANCRW AATEQEEGISLKVECTQLDHEFSCVFAGNPTSCLKLKDERVYWKQVARNLRSQKDICRYSKTAVKTRVCR KDFPESSLKLVSSTLFGNTKPRKEKTEMSPREHIKGKETTPSSLAVTQTMATKAPECVEDPDMANQRKTA LEFCGETWSSLCTFFLSIVQDTSC TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_005130
ORF Size:	702 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_005130.3
RefSeq Size:	1369 bp
RefSeq ORF:	705 bp
Locus ID:	9982
UniProt ID:	Q14512
Cytogenetics:	4p15.32
MW:	26.3 kDa