

## Product datasheet for **RC200820A1V**

### Human **CCDC86 (NM\_024098) AAV Particle**

#### Product data:

Product Type:	AAV Particles
Product Name:	Human CCDC86 (NM_024098) AAV Particle
Tag:	Myc-DDK
Symbol:	CCDC86
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC200820 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGATACACCGTTAAGGCGCAGCCGACGGCTGGGAGGCCAAGGCCGAATCCCCGAGAGCCTCACCT  
CAGTTTCGCGGACGAGACGGGCCCTTGTGGAGTTCGAGTCGAACCCAGAAGAAACGAGGGAGCCCGGGT  
TCCTCCGAGTGTGCAGCGGGCTGGCCTGGGTCCCCGAAAGGCCCGGAAGACAAGCCAGGATCACCC  
CGTCTGCAGCAGGGTGCAGGCTTGGAGTCAACCAAGGCAGCCAGAGCCAGGCGCAGCGTCCCCCAGC  
GTCAGCAAGACCTACACCTGGAGTCGCCTCAAAGACAGCCAGAGTACAGTCTGAATCCCCACGATGTCA  
GCCGAAGCCAAGTGAGGAGGCACCAAAGTGTCTCAGGACCAGGGAGTACTGGCCTCGGAGTTGGCCAG  
AATAAGGAGGAGCTGACCCCGGGGGCCCCCAGCATCAGCTACCGCCGGTCCCAGGATCACCCAGAGCCTT  
ACCCCGGTGAGCAAGCTCCCGGTCCGGAGCCCTCTCAGCCACTACTGGAGCTGACACCCAGGGCACCTGG  
CTCCCCCGGGGTGAGCATGAGCCGAGCAAGCCACCTCCAGCTGGGGAGACGGTGACAGGCGGCTTCGGG  
GCAAAGAAGCGAAAAGTTCTTCATCCCAGGCCCCAGCGTCCAAGAAGTTGAATAAAGAGGAGCTTCCTG  
TAATCCCGAAGGGGAAGCCAAATCGGGCGAGTGTGGAAGGACCGCTCCAAGAAAAGATTCTCCAGAT  
GCTTCAGGACAAGCCCTGCGCACATCGTGGCAGCGGAAGATGAAGGAACGACAGGAGAGGAAGCTGGCC  
AAGGACTTTGCCGTACCTGGAGGAGGAGAAGGAGAGGCGCCGCGCAGGAGAAGAAACAGCCCGGGCTG  
AGAACCTGAAACGCCGCTGGAGAATGAGCGGAAGGCAGAGGTGTCCTCAAGTGATCCGAAACCCCGCAA  
GCTCAAGCGGGCAAAGAAGAAGCAGCTGCGCTCCATTGAGAAGCGGGACACCCTGGCCCTGCTGCAGAAG  
CAGCCGCCCCAGCAGCCGGCAGCCAAGATC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

<b>Protein Sequence:</b>	<p>&gt;RC200820 protein sequence  <span style="color: red;">Red</span>=Cloning site <span style="color: green;">Green</span>=Tags(s)</p> <p>MDTPLRRSRRLGGLRPESPELTSVSRTRRALVEFESNPEETREPGSPPSVQRAGLGSPERPPKTSFGSP  RLQQGAGLESPQGQPEPGAASPQRQDLHLESPQRQPEYSPESPRCQPKPSEEAPKCSQDQGVLAELAQ  NKEELTPGAPQHQLPPVPGSPEPYPGQQAPGPEPSQPLELTPRAPGSPRGQHEPSKPPPAGETVTGGFG  AKKRKGSSSQAPASKKLNKEELPVIPKGGPKSGRVWKRDRSKRFSQMLQDKPLRTSWQRKMKERQERKLA  KDFARHLEEEKERRRQEKKQRRRAENLKRRENERKAEVVQVIRNPAKLRKRAKKQLRSIEKRDTLALLQK  QPPQPAAKI</p> <p><span style="color: red;">TR</span><span style="color: green;">TRPLEQKLI</span><span style="color: green;">SEEDLAANDILDYKDDDDKV</span></p>
<b>Species:</b>	Human
<b>Serotype:</b>	AAV-2
<b>ACCN:</b>	NM_024098
<b>ORF Size:</b>	1080 bp
<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u><a href="#">NM_024098.1</a></u>
<b>RefSeq Size:</b>	2022 bp
<b>RefSeq ORF:</b>	1083 bp
<b>Locus ID:</b>	79080
<b>UniProt ID:</b>	<u><a href="#">Q9H6F5</a></u>
<b>Cytogenetics:</b>	11q12.2
<b>MW:</b>	40.2 kDa