

Product datasheet for **RC200821A1V**

Human CDCA3 (NM_031299) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human CDCA3 (NM_031299) AAV Particle
Tag:	Myc-DDK
Symbol:	CDCA3
Synonyms:	GRCC8; TOME-1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC200821 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGGCTCAGCCAAGAGCGTCCCAGTCACACCAGCGCGGCTCCGCCGACAACAAGCATCTGGCTCGAG
TGGCGGACCCCGTTACCTAGTGCTGGCATCCTGCGCACTCCCATCCAGGTGGAGAGCTCTCCACAGCC
AGGCCTACCAGCAGGGGAGCAACTGGAGGCTCTAAACATGCCAGGACTCAGATCCCCGCTCTCCTACT
CTTGGTATTGCACGGACACCTATGAAGACCAGCAGTGGAGACCCCCAAGCCCACTGGTGAACAGCTGA
GTGAAGTATTTGAACTGAAGACTCTAAATCAAATCTTCCCCAGAGCCTGTTCTGCCCCAGAGGCACC
TTTATCTTCTGAATTGGACTTGCCTCTGGGTACCCAGTTATCTGTTGAGGAACAGATGCCACCTTGGAAC
CAGACTGAGTTCCCTCCAAACAGGTGTTTTCCAAGGAGGAAGCAAGACAGCCACAGAAACCCCTGTGG
CCAGCCAGAGCTCCGACAAGCCCTCAAGGGACCTGAGACTCCCAGATCTTCAGGTTCTATGCGCAATAG
ATGGAACCAACAGCAGCAAGGTACTAGGGAGATCCCCCTCACCATCCTGCAGGATGACAACCTCCCT
GGCACCTGACACTACGACAGGGTAAGCGGCCCTCACCCCTAAGTGAATGTTAGTGAAGTAAAGGAAG
GAGCCATTCTTGAAGTGGACGACTTCTGAAACTGGAGGACGAGCATGGGAGCAAGGCCAGGACCATGA
CAAGGAAAATCAGCACTTCCCTTGGTGGAGAGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC200821 protein sequence Red=Cloning site Green=Tags(s) MGSAKSVPTPARPPPHNKLARVADPRSPSAGILRTPIQVESSPQGPLPAGEQLEGLKHAQSDPRSPT LGIARTPMKTSSGDPPSPLVKQLSEVFETEDSKSNLPPEPVLPEAPLSSELDLPLGTQLSVEEQMPPWN QTEFPSKQVFSKEEARQPTETPVASQSSDKPSRDPETPRSSGSMRNRWKPNSSKVLGRSPLTILQDDNSP GTLTLRQGKRPSPLSENVSELKEGAILGTGRLLKTGGRAWEQGDHDKENQHFLVES TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_031299
ORF Size:	804 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_031299.3</u>
RefSeq Size:	1944 bp
RefSeq ORF:	807 bp
Locus ID:	83461
UniProt ID:	<u>Q99618</u>
Cytogenetics:	12p13.31
MW:	29 kDa