

## Product datasheet for **RC200830L4V**

### **METRN (NM\_024042) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	METRN (NM_024042) Human Tagged ORF Clone Lentiviral Particle
Symbol:	METRN
Synonyms:	C16orf23; c380A1.2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_024042
ORF Size:	879 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200830).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_024042.2</a>
RefSeq Size:	1109 bp
RefSeq ORF:	882 bp
Locus ID:	79006
UniProt ID:	<a href="#">Q9UJH8</a>
Cytogenetics:	16p13.3
Protein Families:	Secreted Protein
MW:	31 kDa



[View online »](#)

**Gene Summary:**

Meteorin regulates glial cell differentiation and promotes the formation of axonal networks during neurogenesis (Nishino et al., 2004 [PubMed 15085178]).[supplied by OMIM, Mar 2008]