

Product datasheet for RC200847L4V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

WNT5B (NM_030775) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Product Name: WNT5B (NM 030775) Human Tagged ORF Clone Lentiviral Particle

Symbol: WNT5B

Mammalian Cell Puromycin

Selection:

Vector: pLenti-C-mGFP-P2A-Puro (PS100093)

Tag: mGFP

ACCN: NM_030775

ORF Size: 1077 bp

ORF Nucleotide

The ORF insert of this clone is exactly the same as(RC200847).

OTI Disclaimer:

Sequence:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeq: <u>NM 030775.2</u>, <u>NP 110402.2</u>

 RefSeq Size:
 2184 bp

 RefSeq ORF:
 1080 bp

 Locus ID:
 81029

 UniProt ID:
 Q9H1J7

 Cytogenetics:
 12p13.33

Domains: wnt

Protein Families: Secreted Protein





WNT5B (NM_030775) Human Tagged ORF Clone Lentiviral Particle - RC200847L4V

Protein Pathways: Basal cell carcinoma, Hedgehog signaling pathway, Melanogenesis, Pathways in cancer, Wnt

signaling pathway

MW: 40.3 kDa

Gene Summary: The WNT gene family consists of structurally related genes which encode secreted signaling

proteins. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. This gene is a member of the WNT gene family. It encodes a protein which shows 94% and 80% amino acid identity to the mouse Wnt5b protein and the human WNT5A protein, respectively.

Alternative splicing of this gene generates 2 transcript variants. [provided by RefSeq, Jul 2008]