

Product datasheet for **RC200856A1V**

Human Protein **FAM195A (MCRIP2) (NM_138418) AAV Particle**

Product data:

Product Type: AAV Particles

Tag: Myc-DDK

Symbol: Protein FAM195A

Synonyms: C16orf14; c349E10.1; FAM195A

Mammalian Cell None

Selection:

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC200856 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGTACACCATCACCAAGGGGCCAGCAAGCTGGTCGCGCAGCGCCGCACAGTCCCACGCAGCAGCAGG
TGGAGGGCCGGCTCGCGAGCTCCTGAAATGCCGCGAGCCCGCGCCGACCTCGCAGCCCCGCGGGC
GCAGCCCTTTGCGCAGCCGCCGGGACCCTGGCCCTGTGAGTCCAGGGCCAAGGCTTGTGTCAATCGT
GTGAATGGCCGGCGGGCCCCCTCCACGTCCCCATCCTTCGAGGGGACCCAGGAGACCTACACAGTGGCCC
ACGAGGAGAATGTCCGCTTTGTGTCCGAAGCCTGGCAGCAGGTGCAACAGCAGCTGGATGGTGGCCAGC
CGGTGAGGGCGGCCAAGGCTGTGCAGTACGTGGAGAGGACCCCAATCCCCGGCTGCAGAACTTTGTG
CCCATTTGACCTAGACGAGTGGTGGCGCAGCAGTTCTGGCGAGAATCACCAGCTGTTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC200856 protein sequence
Red=Cloning site Green=Tags(s)

MYTITKGPSKLVAQRRTGPTQQQVEGRLGELLKCRQPAPPTSQPPRAQPFAPPPGPWPLSSPGPRLVFNR
VNGRRAPSTSPSFEGTQETYTVAHEENVRFVSEAWQQVQQQLDGGPAGEGPRPVQYVERTPNRLQNFV
PIDLDEWQAQFLARITSCS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human

Serotype: AAV-2



ACCN:	NM_138418
ORF Size:	480 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_138418.2</u>
RefSeq Size:	928 bp
RefSeq ORF:	483 bp
Locus ID:	84331
UniProt ID:	<u>Q9BUT9</u>
Cytogenetics:	16p13.3
MW:	17.8 kDa