

Product datasheet for **RC200866A1V**

Human **C11orf73 (HIKESHI) (NM_016401) AAV Particle**

Product data:

Product Type: AAV Particles
Tag: Myc-DDK
Symbol: C11orf73
Synonyms: C11orf73; HLD13; HSPC138; HSPC179; L7RN6; OPI10
Mammalian Cell None
Selection:
Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC200866 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGTTTGGCTGCTTGGTGGCGGGGAGGCTGGTGCAAACAGCTGCACAGCAAGTGGCAGAGGATAAATTTG
TTTTTGACTTACCTGATTATGAAAGTATCAACCATGTTGTGGTTTTATGCTGGGAACAATCCCATTTCC
TGAGGGAATGGGAGGATCTGTCTACTTTTCTTATCCTGATTCAAATGGAATGCCAGTATGGCAACTCCTA
GGATTTGTCACGAATGGGAAGCCAAGTCCATCTTCAAATTTTCAGGTCTTAAATCTGGAGAAGGAAGCC
AACATCCTTTTGGAGCCATGAATATTGTCCGAACCTCATCTGTTGCTCAGATTGGAATTTTCAGTGAATT
ATTAGACAGTATGGCTCAGCAGACTCCTGTAGGTAATGCTGCTGTATCCTCAGTTGACTCATTCACTCAG
TTCACACAAAAGATGTTGGACAATTTCTACAATTTGCTTCATCATTGCTGTCTCTCAGGCCCAGATGA
CACCAAGCCCATCTGAAATGTTTCATTCCGGCAAATGTGGTTCTGAAATGGTATGAAAACCTTTCAAAGACG
ACTAGCACAGAACCCTCTCTTTTGAAAAACA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC200866 protein sequence
Red=Cloning site Green=Tags(s)

MFGCLVAGRLVQTAAQQVAEDKFVFDLPDYESINHVVFVFLGTIPFPEGMGGSVYFSSYPDSNGMPVWQLL
GFVTNGKPSAIFKISGLKSGEGSQHPFGAMNIVRTPSVAQIGISVELLDMAQQTPVGNAAVSSVDSFTQ
FTQKMLDNFYNFASSFAVSQAQMTSPSEMFIPANVVLKWFYENFQRRLAQNPLFWKT

TRTRPLE**Q**KLISEEDLAANDILDYKDDDDKV

Species: Human



Serotype:	AAV-2
ACCN:	NM_016401
ORF Size:	591 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_016401.2
RefSeq Size:	1187 bp
RefSeq ORF:	594 bp
Locus ID:	51501
UniProt ID:	Q53FT3
Cytogenetics:	11q14.2
MW:	21.6 kDa