

Product datasheet for **RC200870A1V**

Human DCUN1D5 (NM_032299) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human DCUN1D5 (NM_032299) AAV Particle
Tag:	Myc-DDK
Symbol:	DCUN1D5
Synonyms:	SCCRO5
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC200870 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCCGGTGAAGAAGAAGAGAAAATCCCCTGGGGTGGCAGCAGCAGTAGCGGAAGACGGAGGCCTCAAAA
AGTGTAATAATCTCCAGCTATTGCAGATCCCAACCCCTGCTAGACTAATAAGTGGAGAGGAACATTTTTC
AAGCAAGAAGTGCCTGGCTTGGTTTTATGAATATGCAGGTCCTGATGAAGTTGTAGGGCCAGAAGGAATG
GAAAAATTTTGTGAAGACATTGGTGTGAACCTGAAAAATTATTATGTTAGTTTTAGCGTGGAATTTGG
AGGCTGAAAGCATGGGATTTTTACCAAGGAAGAATGGTTAAAGGGAATGACTTCATTACAGTGTGACTG
CACAGAAAAGTTACAAAACAAATTTGACTTTTTGCGCTCACAGTTGAATGATATTTTCGTCATTTAAGAAT
ATCTACAGATATGCCTTTGATTTTGCAAGGGATAAAGATCAGAGAAGCCTTGATATTGATCTGCTAAAT
CTATGTTAGCTCTTCTGCTTGGGAGGACATGGCCACTGTTTTTCAGTATTTTACCAGTACCTGGAGCAATC
AAAGTATCGTGTATGAACAAAGATCAATGGTACAATGTATTAGAATTCAGCAGAACAGTCCATGCTGAT
CTTAGTAATATGATGAAGATGGTGCTTGGCCTGTTCTTCTTGATGAATTTGTTGAGTGGCAAAAAGTCC
GTCAGACATCA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC200870 protein sequence Red=Cloning site Green=Tags(s) MPVKKKRKSPGVAAAVAEDGGLKKCKISSYCRSQPPARLISGEEHFSSKKCLAWFYEYAGPDEVVGPEGM EKFCEDIGVEPENIIMLVLAWKLEAESMGFFTKEEWLKGMTSLQCDCTEKLQNKFDLRSQNLNDISSFKN IYRYAFDFARDKDQRLDIDTAKSMLALLLGRTWPLFSVFYQYLEQSKYRVMNKDQWYNVLEFSRTVHAD LSNYDEDGAWPVLLDEFVEWQKVRQTS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_032299
ORF Size:	711 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_032299.2</u>
RefSeq Size:	12732 bp
RefSeq ORF:	714 bp
Locus ID:	84259
UniProt ID:	<u>Q9BTE7</u>
Cytogenetics:	11q22.3
MW:	27.5 kDa