

## Product datasheet for **RC200888L3V**

### **XPNPEP3 (NM\_022098) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	XPNPEP3 (NM_022098) Human Tagged ORF Clone Lentiviral Particle
Symbol:	XPNPEP3
Synonyms:	APP3; ICP55; NPHPL1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_022098
ORF Size:	1521 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200888).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_022098.2</a>
RefSeq Size:	8027 bp
RefSeq ORF:	1524 bp
Locus ID:	63929
UniProt ID:	<a href="#">Q9NQH7</a>
Cytogenetics:	22q13.2
Domains:	Peptidase_M24, AMP_N
Protein Families:	Druggable Genome, Protease



[View online »](#)

**MW:** 57.1 kDa

**Gene Summary:** The protein encoded by this gene belongs to the family of X-pro-aminopeptidases that utilize a metal cofactor, and remove the N-terminal amino acid from peptides with a proline residue in the penultimate position. This protein has been shown to localize to the mitochondria of renal cells, and have a role in ciliary function. Mutations in this gene are associated with nephronophthisis-like nephropathy-1. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene, however, expression of some of these isoforms in vivo is not known.[provided by RefSeq, Mar 2011]