

## Product datasheet for **RC200889L4V**

### NMNAT2 (NM\_170706) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	NMNAT2 (NM_170706) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NMNAT2
Synonyms:	C1orf15; PNAT2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_170706
ORF Size:	906 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200889).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_170706.2</a>
RefSeq Size:	5467 bp
RefSeq ORF:	909 bp
Locus ID:	23057
UniProt ID:	<a href="#">Q9BZQ4</a>
Cytogenetics:	1q25.3
Protein Pathways:	Metabolic pathways, Nicotinate and nicotinamide metabolism
MW:	34 kDa



[View online »](#)

**Gene Summary:**

This gene product belongs to the nicotinamide mononucleotide adenylyltransferase (NMNAT) enzyme family, members of which catalyze an essential step in NAD (NADP) biosynthetic pathway. Unlike the other human family member, which is localized to the nucleus, and is ubiquitously expressed; this enzyme is cytoplasmic, and is predominantly expressed in the brain. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]