

Product datasheet for **RC200890L3V**

DPCD (NM_015448) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | DPCD (NM_015448) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | DPCD |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_015448 |
| ORF Size: | 609 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC200890). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_015448.1 |
| RefSeq Size: | 858 bp |
| RefSeq ORF: | 612 bp |
| Locus ID: | 25911 |
| UniProt ID: | Q9BVM2 |
| Cytogenetics: | 10q24.32 |
| MW: | 23.2 kDa |


[View online »](#)

Gene Summary:

This gene in mouse encodes a protein that may be involved in the generation and maintenance of ciliated cells. In mouse, expression of this gene increases during ciliated cell differentiation, and disruption of this gene has been linked to primary ciliary dyskinesia. [provided by RefSeq, Jul 2016]