

Product datasheet for **RC200899A1V**

Human TRIR (NM_024038) AAV Particle

Product data:

Product Type: AAV Particles
Tag: Myc-DDK
Symbol: TRIR
Synonyms: C19orf43; fSAP18; TERCIR
Mammalian Cell None
Selection:
Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC200899 representing NM_024038
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**C

ATGGCTGCCCGAGGAGACGGGCGGAGCCTCAGGGCCGGGAGGCTCCGGGCCCCGGGGCGGTGGCGGTG
GCGGGAGCCGTTGGGCTGAGTCGGGATCGGGACGTCGCCGAGAGCGGGACGAGGAGGTGTCGGGCGC
GGGTTTCGAGCCCGGTGTCGGGCGGCGTGAACCTTGTTCGCCAACGACGGCAGCTTCCTGGAGCTGTTCAAG
CGGAAGATGGAGGAGGAGCAGCGGCAGCGGAGGAGGCCGCCCGGGTCCGCAGCGACCCGACCACT
CGGCGCCGCCCGCTGGCCCCGGGGATCCGAAGAGGAAGGGCGGTCCGGGCTCCACACTTAGCTTCGTGGG
CAAACGCAGAGGCGGGAACTAGCCCTCAAGACGGGAATAGTAGCCAAGAAGCAGAAGACGGAGGAT
GAGGTATTAACAAGTAAAGGTGACGCGTGGGCCAAGTACATGGCAGAAGTAAAAAGTACAAAGCTCACC
AGTGCGGTGACGATGATAAACTCGGCCCTGGTGAAA

AC**CGGCCGC**TCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGA
TTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC200899 representing NM_024038
Red=Cloning site Green=Tags(s)

MAARGRRAEPQGREAPGPAGGGGGSRWAESGSGTSPESGDEEVSGAGSSPVSGGVNLFANDGSFLELFK
RKMEEEQRQRQEEPPGPQRPDQSAAGPGDPKRKGGPGSTLSFVGKRRGGNKLALKTGIVAKKQKTED
EVLTSKGDWAKYMAEVKKYKAHQCGDDDKTRPLVK

TRRLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2



ACCN:	NM_024038
ORF Size:	528 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_024038.2</u> , <u>NP_076943.1</u>
RefSeq Size:	879 bp
RefSeq ORF:	531 bp
Locus ID:	79002
UniProt ID:	<u>Q9BQ61</u>
Cytogenetics:	19p13.13
MW:	18.4 kDa