

Product datasheet for **RC200900A1V**

Human SSU72 (NM_014188) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human SSU72 (NM_014188) AAV Particle
Tag: Myc-DDK
Symbol: SSU72
Synonyms: HSPC182; PNAS-120
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC200900 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCCGTCGTCCCGCTGCGGGTGGCGGTGGTGTGCTCGAGCAACCAGAACCGGAGCATGGAGGCGCACA
 ACATCCTCAGCAAACGGGGATTACGCGTCCGATCCTTTGGAACAGGGACTCACGTGAAGCTTCCAGGACC
 AGCTCCCGACAAGCCCAATGTTTATGATTTCAAACCACATATGACCAGATGTACAATGATCTTCTTAGG
 AAAGACAAAGAAGCTCTATACACAGAATGGGATTTACATATGCTGGACAGAAATAAGAGAATCAAGCCCC
 GGCCAGAAAGATTCCAGAAGTCAAGACCTGTTTGATCTGATCCTCACTTGCGAAGAGAGAGTGTATGA
 CCAGGTGGTGAAGATCTGAATTCAGAGAACAGGAGACCTGCCAGCCCGTGCAGTGGTCAATGTGGAC
 ATCCAGGACAACCACGAGGAGGCCACCCTGGGGCGTTTCTCATCTGTGAGCTCTGCCAGTGTATCCAGC
 ACACGGAAGACATGGAGAACGAGATCGACGAGCTGCTGCAGGAGTTCGAGGAGAAGAGTGGCCGCACCTT
 TCTGCACACCGTCTGCTTCTAC

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC200900 protein sequence
 Red=Cloning site Green=Tags(s)

MPSSPLRVAVVCSSNQNRSM EAHN ILSKRGFSVRSFGTGTHVKLP GPAPDKPNVYDFKTTYDQMYNDLLR
 KDKELYTQNGILHMLDRNKRIKPRPERFQNC KDLFDLILTCEERVYDQVVEDLNSREQETCQPVHVVNVD
 IQDNHEEATLGAF LICELCQCIQHTEDMENEIDELLQEFEEKSGRTFLHTVCFY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human



[View online »](#)

Serotype:	AAV-2
ACCN:	NM_014188
ORF Size:	582 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_014188.2</u>
RefSeq Size:	1319 bp
RefSeq ORF:	585 bp
Locus ID:	29101
UniProt ID:	<u>Q9NP77</u>
Cytogenetics:	1p36.33
MW:	22.6 kDa