

Product datasheet for **RC200900L1V**

SSU72 (NM_014188) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SSU72 (NM_014188) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SSU72
Synonyms:	HSPC182; PNAS-120
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_014188
ORF Size:	582 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200900).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_014188.2
RefSeq Size:	1319 bp
RefSeq ORF:	585 bp
Locus ID:	29101
UniProt ID:	Q9NP77
Cytogenetics:	1p36.33
Domains:	Ssu72
Protein Families:	Druggable Genome


[View online »](#)

MW: 22.6 kDa

Gene Summary: Protein phosphatase that catalyzes the dephosphorylation of the C-terminal domain of RNA polymerase II. Plays a role in RNA processing and termination. Plays a role in pre-mRNA polyadenylation via its interaction with SYMPK.[UniProtKB/Swiss-Prot Function]