

Product datasheet for **RC200906A1V**

Human SNX11 (NM_152244) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human SNX11 (NM_152244) AAV Particle
Tag:	Myc-DDK
Symbol:	SNX11
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC200906 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGCTTTTGGTGTAGGATGTCGGAGAACCAAGAACAGGAGGAGGTGATTACAGTGCCTGTTTCAGGACC
CCCGAGTGCAGAATGAGGGCTCCTGGAACCTTATGTGGATTATAAGATATTCCTCCATACCAACAGCAA
AGCCTTTACTGCCAAGACTTCCTGTGTGCGGCGCCGCTACCGTGAGTTCGTGTGGCTGAGAAAGCAGCTA
CAGAGAAATGCTGGTTTGGTGCCTGTTCTGAACCTCCTGGGAAGTCAACCTTCTCGGCACCTCAGATG
AGTTCATTGAGAAGCGACGACAAGGTCTGCAGCACTTCCTTGAAAAGGTCTGCAGAGTGTGGTTCTCCT
GTCAGACAGCCAGTTGCACCTATTCCTGCAAAGCCAGCTCTCGGTGCCTGAGATAGAAGCCTGTGTCCAG
GGCCGAAGTACCATGACTGTGTCTGATGCCATTCTTCGATATGCTATGTCAAACCTGTGGCTGGGCCAGG
AAGAGAGGCAGAGCTCTTCTCACCTGGCTAAAGGAGACCAGCCTAAGAGTTGCTGCTTTCTTCCAAGATC
GGGTAGGAGGAGCTCTCCCTCACCGCCTCCCAGTGAAGAAAAGGACCATTTAGAAGTGTGGGCTCCAGTT
GTTGACTCTGAGGTTCTTCTTGAAAAGTCCCACTCTCCACCCCTCTCCTCACCATTATGCTGTGATT
TTGGAAGACCCAAAGAGGGAACCTCCACTCTTCAGTCTGTGAGGAGGGCTGTGGGAGGAGATCATGCTGT
GCCTTTGGACCCTGGTCAGTTAGAAACAGTTTTGAAAAG

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC200906 protein sequence Red=Cloning site Green=Tags(s) MGFWCRMSENQEVEEITVRVQDPRVQNEGSWNSYVDYKIFLHTNSKAFTAKTSCVRRRYREFVWLRKQL QRNAGLVPVPELPGKSTFFGTSTDEFIEKRRQGLQHFLEKVLQSVVLLSDSQLHLFLQSQLSVPEIEACVQ GRSTMTVSDAILRYAMSNCGWAQEERQSSSHLAKGDQPKSCCFLPRSGRRSSPSPPPSEEKDHLEWAPV VDSEVPSLESPTLPLSSPLCCDFGRPKEGTSTLQSVRRAVGGDHAVPLDPGQLETVLEK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_152244
ORF Size:	810 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_152244.1, NP_689450.1</u>
RefSeq Size:	2402 bp
RefSeq ORF:	813 bp
Locus ID:	29916
UniProt ID:	<u>Q9Y5W9</u>
Cytogenetics:	17q21.32
MW:	30.4 kDa