

Product datasheet for **RC200910L4V**

RASL11B (NM_023940) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	RASL11B (NM_023940) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RASL11B
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_023940
ORF Size:	744 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200910).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_023940.2
RefSeq Size:	1979 bp
RefSeq ORF:	747 bp
Locus ID:	65997
UniProt ID:	Q9BPW5
Cytogenetics:	4q12
Domains:	ras, RAS, RHO, RAB
Protein Families:	Stem cell - Pluripotency
MW:	27.5 kDa



[View online »](#)

Gene Summary:

RASL11B is a member of the small GTPase protein family with a high degree of similarity to RAS (see HRAS, MIM 190020) proteins.[supplied by OMIM, Nov 2008]