

Product datasheet for **RC200924L1V**

RAP1GDS1 (NM_021159) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	RAP1GDS1 (NM_021159) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RAP1GDS1
Synonyms:	GDS1; SmgGDS
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_021159
ORF Size:	1821 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200924).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_021159.3 , NP_066982.2
RefSeq Size:	3767 bp
RefSeq ORF:	1824 bp
Locus ID:	5910
UniProt ID:	P52306
Cytogenetics:	4q23
Domains:	Armadillo_seg
MW:	66.3 kDa



[View online »](#)

Gene Summary:

The smg GDP dissociation stimulator (smgGDS) protein is a stimulatory GDP/GTP exchange protein with GTPase activity (Riess et al., 1993 [PubMed 8262526]).[supplied by OMIM, Feb 2010]