

Product datasheet for RC200929A1V

Human IFI35 (NM_005533) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human IFI35 (NM_005533) AAV Particle
Tag: Myc-DDK
Symbol: IFI35
Synonyms: IFP35
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC200929 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGTCAGCCCCACTGGATGCCGCCCTCCACGCCCTTCAGGAGGAGCAGGCCAGACTCAAGATGAGGCTGT
 GGGACCTGCAGCAGCTGAGAAAGGAGCTCGGGGACTCCCCAAAGACAAGGTCCCATTTTCAGTGCCCAA
 GATCCCCCTGGTATTCGAGGACACACCCAGCAGGACCCGGAAGTGCCCTAAGTCTTTAGTTTCCAATTTG
 CGGATCCACTGCCCTCTGCTTGGGGCTCTGCTCTGATCACCTTTGATGACCCCAAAGTGGCTGAGCAGG
 TGCTGCAACAAAAGGAGCACACGATCAACATGGAGGAGTGCCGGCTGCGGGTGCAGGTCCAGCCCTTGA
 GCTGCCCATGGTCACCACTCCAGGTGATGGTGTCCAGCCAGTTGAGTGGCCGGAGGGTGTGGTCACT
 GGATTTCTGCCAGCCTCAGGCTGAGTGAGGAGGAGCTGCTGGACAAGCTAGAGATCTTCTTTGGCAAGA
 CTAGGAACGGAGGTGGCGATGTGGACGTTCCGGAGCTACTGCCAGGGAGTGTCATGCTGGGGTTTGCTAG
 GGATGGAGTGGCTCAGCGTCTGTGCCAAATCGGCCAGTTCACAGTGCCACTGGGTGGGCAGCAAGTCCCT
 CTGAGAGTCTCTCCGTATGTGAATGGGGAGATCCAGAAGGCTGAGATCAGGTGCGAGCCAGTTCCTCGCT
 CGGTACTGGTGTCAACATTCTGATATCTTGGATGGCCGGAGCTGCATGACGTCTGGAGATCCACTT
 CCAGAAGCCCACCCGCGGGGGCGGGAGGTAGAGGCCCTGACAGTCGTACCCCAAGGACAGCAGGGCCTA
 GCAGTCTTCACCTCTGATCAGGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence:	>RC200929 protein sequence Red=Cloning site Green=Tags(s) MSAPLDAALHALQEEQARLKMRLWDLQQLRKELGDSPKDKVPFSVPKIPLVFRGHTQQDPEVPKSLVSNL RIHCPLLAGSALITFDDPKVAEQVLQQKEHTINMEECRLRVQVQPLELPMVTTIQVMVSSQLSGRRVLVT GFPASLRLSEEELLDKLEIFFGKTRNGGGDQVDVRELLPGSVMLGFARDGVAQRLCQIGQFTVPLGGQQVP LRVSPYVNGEIQKAEIRSQPVPRSVLVNIPDILDGPELHDVLEIHFQKPTRGGGEVEALTVVPQGGQGL AVFTSESG TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_005533
ORF Size:	864 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_005533.2 , NP_005524.1
RefSeq Size:	1260 bp
RefSeq ORF:	867 bp
Locus ID:	3430
UniProt ID:	P80217
Cytogenetics:	17q21.31
MW:	31.7 kDa