

## Product datasheet for **RC200930L2V**

### ASMTL (NM\_004192) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ASMTL (NM_004192) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ASMTL
Synonyms:	ASMTLX; ASMTLY; ASTML
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_004192
ORF Size:	1863 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200930).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_004192.2</a>
RefSeq Size:	2149 bp
RefSeq ORF:	1866 bp
Locus ID:	8623
UniProt ID:	<a href="#">O95671</a>
Cytogenetics:	X;Y
Domains:	Methyltransf_2, Maf
MW:	68.8 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene has an N-terminus that is similar to the multicopy associated filamentation (maf) protein of *Bacillus subtilis* and to orfE of *Escherichia coli*, while the C-terminus is similar to N-acetylserotonin O-methyltransferase. This gene is located in the pseudoautosomal region 1 (PAR1) of X and Y chromosomes. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Mar 2010]