

Product datasheet for **RC200932A1V**

Human **CCDC92 (NM_025140)** AAV Particle

Product data:

Product Type: AAV Particles
Tag: Myc-DDK
Symbol: CCDC92
Mammalian Cell None
Selection:
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC200932 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGACTTCACCACATTTCTCGAGTTACGATGAAGGTCCTCTGGATGTCAGCATGGCAGCCACAACTGG
AGAACCAGCTGCACAGCGCACAGAAGAACCTCCTGTTCTTCAGCGGGAGCATGCCAGCAGCTCAAGGG
GCTGCACTCCGAGATCAGGCGGCTGCAGCAGCACTGCACAGATTTAACATATGAGCTGACAGTCAAAAGT
TCGGAACAGACAGGAGACGGGACTTCTAAAGCAGTGAATTAAGAAAAGATGTGAAGAGCTGGAAGCCC
AACTGAAAGTGAAAGAGAACGAAAATGCTGAGTTGTTGAAAGAACTGGAGCAGAAAAACGCGATGATCAC
AGTGCTGGAGAACACCATCAAGGAGCGAGAGAAGAAAGTACCTGGAGGAGCTGAAGGCCAAGAGTCACAAG
CTGACCCTGCTGTCTAGCGAGCTGGAGCAGCGGGCCAGCACCATCGCCTACCTGACCTCCAGCTGCACG
CCGCCAAGAAGAAGCTCATGAGCTCCAGCGGGACCTCAGATGCCAGCCCGTCAGGGAGCCCCGTGCTGGC
CAGCTACAAGCCAGCGCCCCCAAAGACAAGCTACCCGAAACGCCTCGCCGCCGCATGAAAAAGAGCCTC
TCAGCCCCCTTGACCCGGAATTTGAAGAGGTCTACAGATTCGGGGCAGAGAGCAGGAACTCCTTTTGC
GGGAACCAGTGGATGCTATGCCGACCCACCCATTCTGCTGGCTAGGGAGTCCGCCGAGGTCCACCT
CATCAAAGAGAGGCCCTCGTCATCCCCCATCGCCTCCGACCGAAGCGGCGAGCAGCACAGCCGGCC
CGCGAAAAGCCGCACAAGGCCACGTCGGGGTGGCACATCGGATCCACCACGCCACCCGCGCAGGCC
AGCCCCGAGGTGAAGACCTGGCGGTGACCAAGTGAACGGAGGCAAGGTGGTGAAGAACTCAGGGAC
GGACAGAACTGTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATACGACGATAAGGTTAA



Protein Sequence:	>RC200932 protein sequence Red=Cloning site Green=Tags(s) MTSPHFSSYDEGPLDVSMATNLENQLHSAQKNLLFLQREHASTLKGLHSEIRRLQQHCTDLTYELTVKS SEQTGDGTSKSSSELKKRCEELEAQLKVENENAELLKELEQKNAMITVLENTIKEREKKYLEELKAKSHK LTLLSSELEQRASTIAYLTSQLHAAKKKLMSSSGTSDASPSGSPVLASYKPAPPKDKLPETPRRRMKKSL SAPLHPEFEEVYRFGAESRKLLLREPVDAMPDPTPFLLAESAIEVHLIKERPLVIPPIASDRSGEQHSPA REKPHKAHVGVAHRIHHATPPQAQPEVKTLAVDQVNGGKVVRKHSGTDRTV TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_025140
ORF Size:	993 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_025140.1
RefSeq Size:	1814 bp
RefSeq ORF:	996 bp
Locus ID:	80212
UniProt ID:	Q53HC0
Cytogenetics:	12q24.31
MW:	37 kDa