

Product datasheet for **RC200936L1V**

NUDT21 (NM_007006) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | NUDT21 (NM_007006) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | NUDT21 |
| Synonyms: | CFIM25; CPSF5 |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-Myc-DDK (PS100064) |
| Tag: | Myc-DDK |
| ACCN: | NM_007006 |
| ORF Size: | 681 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC200936). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_007006.2 |
| RefSeq Size: | 4407 bp |
| RefSeq ORF: | 684 bp |
| Locus ID: | 11051 |
| UniProt ID: | O43809 |
| Cytogenetics: | 16q13 |
| Domains: | NUDIX |
| MW: | 26.2 kDa |



[View online »](#)

Gene Summary:

The protein encoded by this gene is one subunit of a cleavage factor required for 3' RNA cleavage and polyadenylation processing. The interaction of the protein with the RNA is one of the earliest steps in the assembly of the 3' end processing complex and facilitates the recruitment of other processing factors. This gene encodes the 25kD subunit of the protein complex, which is composed of four polypeptides. [provided by RefSeq, Jul 2008]