

Product datasheet for **RC200941L4V**

FKRP (NM_024301) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | FKRP (NM_024301) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | FKRP |
| Synonyms: | FKTR; LGMD2I; LGMDR9; MDC1C; MDDGA5; MDDGB5; MDDGC5 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_024301 |
| ORF Size: | 1485 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC200941). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_024301.3 |
| RefSeq Size: | 3422 bp |
| RefSeq ORF: | 1488 bp |
| Locus ID: | 79147 |
| UniProt ID: | Q9H9S5 |
| Cytogenetics: | 19q13.32 |
| Protein Families: | Transmembrane |
| MW: | 54.6 kDa |



[View online »](#)

Gene Summary:

This gene encodes a protein which is targeted to the medial Golgi apparatus and is necessary for posttranslational modification of dystroglycan. Mutations in this gene have been associated with congenital muscular dystrophy, cognitive disability, and cerebellar cysts. Several alternatively spliced transcript variants of this gene have been described, but the full-length nature of some of these variants has not been determined. [provided by RefSeq, Oct 2008]