

Product datasheet for **RC200946L4V**

ISG20L2 (NM_030980) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ISG20L2 (NM_030980) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ISG20L2
Synonyms:	FLJ12671
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_030980
ORF Size:	1059 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200946).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_030980.1 , NP_112242.1
RefSeq Size:	2051 bp
RefSeq ORF:	1062 bp
Locus ID:	81875
Cytogenetics:	1q23.1
Domains:	EXOIII
MW:	39.2 kDa


[View online »](#)

Gene Summary:

This gene encodes a 3'-5' exoribonuclease that may be involved in the processing of the 12S pre-rRNA. Pseudogenes have been identified on chromosomes 6 and 11. [provided by RefSeq, Dec 2014]