

## Product datasheet for **RC200952A1V**

### Human NUDT9 (NM\_198038) AAV Particle

#### Product data:

|                |                            |
|----------------|----------------------------|
| Product Type:  | AAV Particles              |
| Tag:           | Myc-DDK                    |
| Symbol:        | NUDT9                      |
| Synonyms:      | NUDT10                     |
| Mammalian Cell | None                       |
| Selection:     |                            |
| Vector:        | pAAV-AC-Myc-DDK (PS100089) |

**ORF Nucleotide Sequence:** >RC200952 ORF sequence  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGCATCGCC**

ATGGCGGGACGCCTCCTGGGAAAGGCTTTAGCCGCGGTGTCTCTCTCTGTCCTTGGCCTTGGCCTCTGTGACTA  
TCAGGTCCTCGCGCTGCCGCGGCATCCAGGCGTTCAGAACTCGTTTTTCATCTTCTTGGTTTCATCTTAA  
TACCAACGTATGTCTGGTTCTAATGGTTCCAAAGAAAATTCTCACAATAAGGCTCGGACGTCTCCTTAC  
CCAGGTTCAAAAGTTGAACGAAGCCAGGTTCTAATGAGAAAGTGGGCTGGCTTGTGAGTGGAAGACT  
ATAAGCCTGTGGAATACACTGCAGTCTCTGTCTGGCTGGACCCAGGTGGGCAGATCCTCAGATCAGTGA  
AAGTAATTTTCTCCCAAGTTTAACGAAAAGGATGGGCATGTTGAGAGAAAGAGCAAGAATGGCCTGTAT  
GAGATTGAAAATGGAAGACCGAGAAATCCTGCAGGACGGAAGTGGTGGGCCGGGGCTTTTGGGGC  
GATGGGGCCCAATCACGCTGCAGATCCCATTATAACCAGATGGAAGGGATAGCAGTGGAATAAAAT  
CATGCATCCTGTTCTGGGAAGCATATCTTACAATTTGTTGCAATAAAAGGAAAGACTGTGGAGAATGG  
GCAATCCAGGGGGGATGGTGGATCCAGGAGAGAAGATTAGTGCCACACTGAAAAGAGAATTTGGTGAGG  
AAGCTCTCAACTCCTTACAGAAAACAGTGCTGAGAAGAGAGAAAATAGAGGAAAAGTTGCACAACTCTT  
CAGCCAAGACCACCTAGTGATATATAAGGGATATGTTGATGATCCTCGAAACACTGATAATGCATGGATG  
GAGACAGAAGCTGTGAATACCATGACGAAACAGGTGAGATAATGGATAATCTTATGCTAGAAGCTGGAG  
ATGATGCTGGAAGAGTGAATGGGTGGACATCAATGATAAACTGAAGCTTTATGCCAGTCACTCTCAATT  
CATCAAACCTTGTGGCTGAGAAACGAGATGCACACTGGAGCGAGGACTCTGAAGCTGACTGCCATGCGTTG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Protein Sequence:</b> | <p>&gt;RC200952 protein sequence</p> <p>Red=Cloning site Green=Tags(s)</p> <p>MAGRLLGKALAAVSLSLALASVTIRSSRCRGIQAFRNSFSSSWFHLNTNVMMSGNSKENSHNKARTSPY<br/> PGSKVERSQVPNEKVGWLVEWQDYKPVEYTAVSVLAGPRWADPQISESNFSPKFNEKDGHVERKSKNGLY<br/> EIENGRPRNPAGRTGLVGRLLGRWGNHAADPIITRWKRDSSGNKIMHPVSGKHILQFVAIKRKDCGEW<br/> AIPGGMVDPGEKISATLKREFGEEALNSLQKTSAEKREIEEKLHKLFSDHLVIYKGYVDDPRNTDNAWM<br/> ETEAVNYHDETGEIMDNLMEAGDDAGVKWVDINDKLKLYASHSQFIKLVAEKRDAHWSEADCHAL</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p> |
| <b>Species:</b>          | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Serotype:</b>         | AAV-2                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>ACCN:</b>             | NM_198038                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>ORF Size:</b>         | 1050 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Buffer:</b>           | PBS with 0.001% Pluronic F68                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Stability:</b>        | AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.                            |
| <b>RefSeq:</b>           | <u><a href="#">NM_198038.1</a></u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>RefSeq Size:</b>      | 2458 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>RefSeq ORF:</b>       | 903 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Locus ID:</b>         | 53343                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>UniProt ID:</b>       | <u><a href="#">Q9BW91</a></u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Cytogenetics:</b>     | 4q22.1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>MW:</b>               | 39.1 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |