

Product datasheet for **RC200964L3V**

Calpain 1 (CAPN1) (NM_005186) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Calpain 1 (CAPN1) (NM_005186) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Calpain 1
Synonyms:	CANP; CANP1; CANPL1; muCANP; muCL; SPG76
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_005186
ORF Size:	2142 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200964).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_005186.2
RefSeq Size:	3110 bp
RefSeq ORF:	2145 bp
Locus ID:	823
UniProt ID:	P07384
Cytogenetics:	11q13.1
Domains:	Calpain_III, EFh
Protein Families:	Druggable Genome, Protease



[View online »](#)

Protein Pathways: Alzheimer's disease, Apoptosis

MW: 81.9 kDa

Gene Summary: The calpains, calcium-activated neutral proteases, are nonlysosomal, intracellular cysteine proteases. The mammalian calpains include ubiquitous, stomach-specific, and muscle-specific proteins. The ubiquitous enzymes consist of heterodimers with distinct large, catalytic subunits associated with a common small, regulatory subunit. This gene encodes the large subunit of the ubiquitous enzyme, calpain 1. Several transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Nov 2010]